


Laboratory Underslung Sink Bowls

Kinetic offer a range of laboratory underslung stainless steel sink bowls that are designed to be easily fitted to all types of laboratory worktop. The range is available in a choice of sizes, shapes, colours and either 304 or 316 quality stainless steel.




Underslung bowl with fluted drainer on sink unit




Underslung bowl dimensions




Underslung bowl without fluted drainer on sink unit



Underslung bowl with fluted drainer on sink unit



Underslung bowl without fluted drainer on sink unit



Underslung bowl with fluted drainer on sink unit



Underslung bowl without fluted drainer on sink unit



Underslung bowl with fluted drainer on sink unit



Underslung bowl without fluted drainer on sink unit

Sizes

Code	Length	Width	Depth
	mm	mm	mm
FPB0035	356	229	180
FPB0036	356	305	160
FPB0037	432	356	180
FPB0038	508	356	178
FPB0038D	508	356	250
FPB0040	533	381	200
FPB0040D	533	381	250
FPB0042	610	457	305

Colour Range

Stainless Steel



Stainless
Steel

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.