




## Laboratory Inset Sink & Drainers

Kinetic offer a range of laboratory inset stainless steel sink & drainers that are designed to be easily fitted to all types of laboratory worktop. Complete with sloping drainer that can be right handed, left handed or on both sides. The range is available in a choice of sizes, shapes, colours and either 304 or 316 quality stainless steel.



Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer dimensions




Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer on sink unit



Single bowl with left hand drainer dimensions




Single bowl with left hand drainer on sink unit




Single bowl with left hand drainer on sink unit



Single bowl with left hand drainer on sink unit




Single bowl with left hand drainer on sink unit



Single bowl with double drainer dimensions




Single bowl double drainer on sink unit



Single bowl double drainer on sink unit



Single bowl double drainer on sink unit



Single bowl double drainer on sink unit

## Sizes

| Code    | Length | Width | Depth | Bowl Length | Bowl Width | Drainer    |
|---------|--------|-------|-------|-------------|------------|------------|
|         | mm     | mm    | mm    | mm          | mm         | Drainer    |
| FPI0012 | 965    | 525   | 178   | 432         | 356        | Right Hand |
| FPI0013 | 965    | 525   | 178   | 432         | 356        | Left Hand  |
| FPI0014 | 1030   | 500   | 250   | 508         | 356        | Right Hand |
| FPI0015 | 1030   | 500   | 250   | 508         | 356        | Right Hand |
| FPI0016 | 1200   | 650   | 300   | 610         | 457        | Left Hand  |
| FPI0017 | 1200   | 650   | 300   | 610         | 457        | Right Hand |
| FPI0018 | 1376   | 525   | 178   | 432         | 356        | Double     |
| FPI0019 | 1448   | 525   | 250   | 508         | 356        | Double     |

## Colour Range

### Stainless Steel



Stainless Steel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.