




## Laboratory Inset Sink & Drainers


Kinetic offer a range of laboratory inset stainless steel sink & drainers that are designed to be easily fitted to all types of laboratory worktop. Complete with sloping drainer that can be right handed, left handed or on both sides. The range is available in a choice of sizes, shapes, colours and either 304 or 316 quality stainless steel.



Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer dimensions




Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer on sink unit



Single bowl with right hand drainer on sink unit



Single bowl with left hand drainer dimensions




Single bowl with left hand drainer on sink unit




Single bowl with left hand drainer on sink unit



Single bowl with left hand drainer on sink unit




Single bowl with left hand drainer on sink unit



Single bowl with double drainer dimensions




Single bowl double drainer on sink unit



Single bowl double drainer on sink unit



Single bowl double drainer on sink unit



Single bowl double drainer on sink unit

## Sizes

Code	Length	Width	Depth	Bowl Length	Bowl Width	Drainer
	mm	mm	mm	mm	mm	Drainer
FPI0012	965	525	178	432	356	Right Hand
FPI0013	965	525	178	432	356	Left Hand
FPI0014	1030	500	250	508	356	Right Hand
FPI0015	1030	500	250	508	356	Right Hand
FPI0016	1200	650	300	610	457	Left Hand
FPI0017	1200	650	300	610	457	Right Hand
FPI0018	1376	525	178	432	356	Double
FPI0019	1448	525	250	508	356	Double

## Colour Range

### Stainless Steel



Stainless Steel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.