

Laboratory Free Standing Table

Kinetic laboratory free-standing tables and benches are a flexible independent laboratory table that compliments any scientific laboratory. Used in Biology, chemistry and physics classrooms and laboratories across the industry, Kinetic laboratory free-standing tables and benches come in a variety of worktop materials and a bespoke range of sizes and colours.

Constructed with a fully welded steel frame, adjustable feet for levelling and a choice of central or back bottom rail. The laboratory free-standing tables and benches can take a static load as standard of up to 125kg, with the option of heavy-duty feet for a greater load bearing.

The steel frame is epoxy powder coated and available in a range of colours.



Free Standing




Free Standing



Free Standing



Free Standing



Free Standing



Free Standing



Free Standing

Dimensions

Length	600 - 1200mm
Depth	600 - 900mm
Height	750 - 950mm

Colour Range


Laboratory Worktop - Trespa Toplab® Base



Laboratory Worktop - Trespa Toplab® Plus



Laboratory Worktop - Laminate



Steelwork



Colours on this page may not be a true representation.
For accurate colour samples, please contact us.