

## Laboratory Mobile Table

Kinetic laboratory mobile tables and benches are a flexible independent laboratory table that compliments any scientific laboratory. Used in Biology, chemistry and physics classrooms and laboratories across the industry, Kinetic laboratory mobile tables and benches come in a variety of worktop materials and a bespoke range of sizes and colours.

Constructed with a fully welded steel frame, swivel castors and a choice of central or back bottom rail. The laboratory mobile tables and benches can take a static load as standard of up to 125kg, with the option of using heavy-duty Castors for a greater load bearing.


The steel frame is epoxy powder coated and available in a range of colours.



Mobile Table



Mobile Table



Mobile Table



Mobile Table



Mobile Table



Mobile Table

## Dimensions

|        |              |
|--------|--------------|
| Length | 600 - 1200mm |
| Depth  | 600 - 900mm  |
| Height | 750 - 950mm  |

## Colour Range


### Laboratory Worktop - Trespa Toplab® Base



### Laboratory Worktop - Trespa Toplab® Plus



### Laboratory Worktop - Laminate



### Steelwork



Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.