

## Modular Mobile Tables & benches

Kinetic Laboratories modular mobile tables and benches are a flexible independent laboratory table, available in a standard size of 1500mm. The modular mobile laboratory table is designed to be fully compatible with Kinetic modular suspended units. The scientific modular mobile table comes complete with a Trespa ® Toplab ® BASE available in a variety of colours.


Constructed with a fully welded steel frame and swivel castors. The modular mobile tables and benches can take a static load as standard of up to 125kg, with the option of using heavy-duty castors for a greater load bearing.

The steel frame is epoxy powder coated and available in a range of colours.

Each table is able to support up to 3 modular storage units. The modular storage units can be purchased at the same time or they can be easily installed into an existing Kinetic Laboratories modular mobile table at a later date.



Modular scientific  
Mobile table



Modular scientific  
Mobile table



Modular scientific  
Mobile table



Modular scientific  
Mobile table



Modular scientific  
Mobile table



Modular scientific  
Mobile table

## Associated Products



Mobile table with  
units



Mobile table with  
unit



Mobile table with  
unit



Mobile table with  
units



Mobile table with  
units























Mobile table with  
units

## Dimensions

Length (mm)	1500
Depth (mm)	750
Height (mm)	900

## Colour Range

Laboratory Worktop - Trespa Toplab® Base											
	White		Sand		Icy Blue		Pastel Grey		Silver Grey		Mid Grey
	Speckle White		Speckle Pastel Grey		Speckle Silver Grey		Speckle Icy Blue		Speckle Powder Blue		Speckle Aquamarine
	Speckle Sand		Oak								
Steelwork											
	Ral 9010 White		Ral 7035 Light Grey		Ral 7012 Basalt Grey		Ral 7021 Dark Grey		Ral 5003 Sapphire Blue		Ral 5005 Signal Blue

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.