

Modular Mobile Tables & benches

Kinetic Laboratories modular mobile tables and benches are a flexible independent laboratory table, available in a standard size of 1500mm. The modular mobile laboratory table is designed to be fully compatible with Kinetic modular suspended units. The scientific modular mobile table comes complete with a Trespa® Toplab® BASE available in a variety of colours.

Constructed with a fully welded steel frame and swivel castors. The modular mobile tables and benches can take a static load as standard of up to 125kg, with the option of using heavy-duty castors for a greater load bearing.

The steel frame is epoxy powder coated and available in a range of colours.


Each table is able to support up to 3 modular storage units. The modular storage units can be purchased at the same time or they can be easily installed into an existing Kinetic Laboratories modular mobile table at a later date.



Modular scientific Mobile table



Modular scientific Mobile table



Modular scientific Mobile table



Modular scientific Mobile table



Modular scientific Mobile table



Modular scientific Mobile table

Associated Products



Mobile table with units



Mobile table with unit



Mobile table with unit



Mobile table with units



Mobile table with units



Mobile table with units

Dimensions

Length (mm)	1500
Depth (mm)	750
Height (mm)	900

Colour Range

Laboratory Worktop - Trespa Toplab® Base



Steelwork



Colours on this page may not be a true representation.
For accurate colour samples, please contact us.