

## Modular Free Standing Tables & benches

Kinetic Laboratories modular free-standing tables and benches are a flexible independent laboratory table, available in a standard size of 1500mm. The modular free-standing laboratory table is designed to be fully compatible with Kinetic modular suspended units. The scientific modular free-standing table comes complete with a Trespa ® Toplab ® BASE available in a variety of colours.

Constructed with a fully welded steel frame and swivel castors. The modular free-standing tables and benches can take a static load as standard of up to 125kg, with the option of using heavy-duty feet for a greater load bearing.

The steel frame is epoxy powder coated and available in a range of colours.

Each table is able to support up to 3 modular storage units. The modular storage units can be purchased at the same time or they can be easily installed into an existing Kinetic Laboratories modular mobile table at a later date.



Modular scientific free standing table



Modular scientific free standing table




Modular scientific free standing table



Modular scientific free standing table



Modular scientific free standing table



Modular scientific free standing table

## Associated Products



Free Standing table with unit



Free Standing table with units



Free Standing table with unit



Free Standing table with Units



Free Standing table with units



Free Standing table with units

## Dimensions

|             |      |
|-------------|------|
| Length (mm) | 1500 |
| Depth (mm)  | 750  |
| Height (mm) | 900  |

## Colour Range

### Laboratory Shelf - Trespa Toplab® Base

|  |               |   |                     |   |                     |   |                  |   |                     |   |                    |
|--|---------------|---|---------------------|---|---------------------|---|------------------|---|---------------------|---|--------------------|
|  | White         |  | Sand                |  | Icy Blue            |  | Pastel Grey      |  | Silver Grey         |  | Mid Grey           |
|  | Speckle White |  | Speckle Pastel Grey |  | Speckle Silver Grey |  | Speckle Icy Blue |  | Speckle Powder Blue |  | Speckle Aquamarine |
|  | Speckle Sand  |  | Beech               |  | Oak                 |   |                  |   |                     |   |                    |

### Steelwork

|  |                |   |                     |   |                      |   |                    |   |                        |   |                      |
|--|----------------|---|---------------------|---|----------------------|---|--------------------|---|------------------------|---|----------------------|
|  | Ral 9010 White |  | Ral 7035 Light Grey |  | Ral 7012 Basalt Grey |  | Ral 7021 Dark Grey |  | Ral 5003 Sapphire Blue |  | Ral 5005 Signal Blue |
|--|----------------|---|---------------------|---|----------------------|---|--------------------|---|------------------------|---|----------------------|

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.