

## Modular Free Standing Tables & benches

Kinetic Laboratories modular free-standing tables and benches are a flexible independent laboratory table, available in a standard size of 1500mm. The modular free-standing laboratory table is designed to be fully compatible with Kinetic modular suspended units. The scientific modular free-standing table comes complete with a Trespa ® Toplab ® BASE available in a variety of colours.

Constructed with a fully welded steel frame and swivel castors. The modular free-standing tables and benches can take a static load as standard of up to 125kg, with the option of using heavy-duty feet for a greater load bearing.

The steel frame is epoxy powder coated and available in a range of colours.

Each table is able to support up to 3 modular storage units. The modular storage units can be purchased at the same time or they can be easily installed into an existing Kinetic Laboratories modular mobile table at a later date.



Modular scientific free standing table



Modular scientific free standing table




Modular scientific free standing table



Modular scientific free standing table



Modular scientific free standing table



Modular scientific free standing table

## Associated Products



Free Standing table with unit



Free Standing table with units



Free Standing table with unit



Free Standing table with Units



Free Standing table with units



Free Standing table with units

## Dimensions

Length (mm)	1500
Depth (mm)	750
Height (mm)	900

## Colour Range

### Laboratory Shelf - Trespa Toplab® Base

	White		Sand		Icy Blue		Pastel Grey		Silver Grey		Mid Grey
	Speckle White		Speckle Pastel Grey		Speckle Silver Grey		Speckle Icy Blue		Speckle Powder Blue		Speckle Aquamarine
	Speckle Sand		Beech		Oak						

### Steelwork

	Ral 9010 White		Ral 7035 Light Grey		Ral 7012 Basalt Grey		Ral 7021 Dark Grey		Ral 5003 Sapphire Blue		Ral 5005 Signal Blue
--	----------------	---	---------------------	---	----------------------	---	--------------------	---	------------------------	---	----------------------

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.