



Modular Suspended Storage Units

Kinetic laboratory suspended units can be installed onto the appropriate laboratory steelwork with the ability of being removed if purpose alters or for access to rear services or cleaning. The laboratory furniture range is available in a bespoke choice of sizes and colours.



Single 1 drawer 1 door unit



Single Cupboard unit



Single 1 drawer 1 door unit



Single 1 drawer 1 door unit



Single Cupboard unit



Single Cupboard unit



Single 4 drawer unit



Single 4 drawer unit



Single Cupboard unit



Single 1 drawer 1 door unit



Single 4 drawer unit



Single 4 drawer unit

Associated Products



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table

Dimensions

Overall Height	698mm
Length	450mm
Carcass Depth	500mm

Hardware

Hinges	165° opening hinge
Drawer Sides	White epoxy coated metal drawer sides

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDf)

	White K101PE		Grey U112PE
--	------------------------	---	-----------------------

Laboratory Fronts - Laminate Faced MDF (LFMDf)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.