



## Modular Suspended Storage Units

Kinetic laboratory suspended units can be installed onto the appropriate laboratory steelwork with the ability of being removed if purpose alters or for access to rear services or cleaning. The laboratory furniture range is available in a bespoke choice of sizes and colours.



Single 1 drawer 1 door unit



Single Cupboard unit



Single 1 drawer 1 door unit



Single 1 drawer 1 door unit



Single Cupboard unit



Single Cupboard unit



Single 4 drawer unit



Single 4 drawer unit



Single Cupboard unit



Single 1 drawer 1 door unit



Single 4 drawer unit



Single 4 drawer unit

## Associated Products



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Mobile Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table



Units mounted in Free Standing Table

## Dimensions

Overall Height	698mm
Length	450mm
Carcass Depth	500mm

## Hardware














Hinges	165° opening hinge
Drawer Sides	White epoxy coated metal drawer sides

## Colour Range

### Laboratory Carcass - Melamine Faced MDF (MFMDf)

 <b>White</b> K101PE	 <b>Grey</b> U112PE
---	---

### Laboratory Fronts - Laminate Faced MDF (LFMDf)

 <b>Aluminium</b> D1309WF	 <b>Stone Oak</b> D2217PR	 <b>Natural Maple</b> D375PR	 <b>Beech</b> D381PR	 <b>White</b> K101PE	 <b>Grey</b> U112PE
 <b>Blue</b> U125PE	 <b>Yellow</b> U131PE	 <b>Orange</b> U136PE	 <b>Coral Red</b> U148PE	 <b>Anthracite</b> U164PE	 <b>Black</b> U190PE
 <b>Lime</b> U2504PE	 <b>Caribbean Blue</b> U2510PE	 <b>Magnolia</b> U520PE			

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.