

A Frame

Kinetic A Frame laboratory support system allows independent laboratory worktop support without the use of laboratory under bench storage.

The span frames and the low level back rails can vary in length and will form a series of individual spaces. Allowing the client the option of adding laboratory suspended units, laboratory free standing units or laboratory mobile units.

The storage spaces can also be customised to take clients own equipment such as fridges and freezers or can be left open to provide knee spaces for write-up areas.

Steel frame are epoxy powder coated and are available in a range of colours.



A Frame laboratory support system with Trespa Toplab base worktop



A Frame laboratory support system with Trespa Toplab base worktop



A Frame laboratory support system with Trespa Toplab base worktop



A Frame laboratory support system with Trespa Toplab base worktop



A Frame laboratory support system with Trespa Toplab base worktop



A Frame laboratory support system with Trespa Toplab base worktop

Associated Products



Laboratory mobile units in A frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in A frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in A frame laboratory support system with Trespa Toplab base worktop



Laboratory suspended units in A frame laboratory support system with Trespa Toplab base worktop



Laboratory suspended units in A frame laboratory support system with Trespa Toplab base worktop



Laboratory suspended units in A frame laboratory support system with Trespa Toplab base worktop

Colour Range

Laboratory Worktop - Trespa Toplab® Base

	White		Sand		Icy Blue		Pastel Grey		Silver Grey		Mid Grey
	Speckle White		Speckle Pastel Grey		Speckle Silver Grey		Speckle Icy Blue		Speckle Powder Blue		Speckle Aquamarine
	Speckle Sand		Oak								

Laboratory Worktop - Trespa Toplab® Plus

	White		Pastel Grey		Silver Grey		Slate Grey		Glacier Blue		Pastel Grey / Silver
	Black										

Laboratory Worktop - Epoxy Resin

	Grey		White		Black
--	------	---	-------	---	-------

Laboratory Worktop - Velstone

	Soft White		Silver Ash		Packice		Silver Dust
--	------------	---	------------	---	---------	---	-------------

Steelwork

	Ral 9010 White		Ral 7035 Light Grey		Ral 7012 Basalt Grey		Ral 7021 Dark Grey		Ral 5003 Sapphire Blue		Ral 5005 Signal Blue
--	----------------	---	---------------------	---	----------------------	---	--------------------	---	------------------------	---	----------------------

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.