



Pedestal

Kinetic Pedestal laboratory support system is fixed into position and directly supports the laboratory worktop.

The laboratory under bench unit stands on a plinth in a fixed position and can vary in length.

The storage spaces can be customised using a combination of drawers, shelves and cupboard space.

The Pedestal laboratory support system is available in a range of colours.



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop



Laboratory pedestal units with Trespa Toplab base worktop

Available On


- Melamine Faced MDF
- Laminate MDF

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDf)



Laboratory Fronts - Laminate Faced MDF (LFMDf)



Laboratory Worktop - Trespa Toplab® Base




Laboratory Worktop - Trespa Toplab® Plus



Laboratory Worktop - Epoxy Resin



Laboratory Worktop - Velstone



Colours on this page may not be a true representation.
For accurate colour samples, please contact us.