

## T Frame

Kinetic T Frame laboratory support system allows independent laboratory worktop support with unrestricted continuous space for the user.

The span frames and low level back rails can vary in length and will form a continual unrestricted space. Allowing the client the option of adding laboratory suspended units, laboratory free standing units or laboratory mobile units.

The storage space can also be customised to take clients own laboratory equipment such as fridges and freezers or can be left open to provide knee spaces for write-up areas.


The laboratory steel frames are epoxy powder coated and available in a range of colours.



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop



T Frame laboratory support system with Trespa Toplab base worktop

# Associated Products



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop



Laboratory mobile units in T frame laboratory support system with Trespa Toplab base worktop

# Colour Range

## Laboratory Worktop - Trespa Toplab® Base

	White		Sand		Icy Blue		Pastel Grey		Silver Grey		Mid Grey
	Speckle White		Speckle Pastel Grey		Speckle Silver Grey		Speckle Icy Blue		Speckle Powder Blue		Speckle Aquamarine
	Speckle Sand		Oak								

## Laboratory Worktop - Trespa Toplab® Plus

	White		Pastel Grey		Silver Grey		Slate Grey		Glacier Blue		Pastel Grey / Silver
	Black										

## Laboratory Worktop - Epoxy Resin

	Grey		White		Black
--	------	---	-------	---	-------

## Laboratory Worktop - Velstone

	Soft White		Silver Ash		Packice		Silver Dust
--	------------	---	------------	---	---------	---	-------------

## Steelwork

	Ral 9010 White		Ral 7035 Light Grey		Ral 7012 Basalt Grey		Ral 7021 Dark Grey		Ral 5003 Sapphire Blue		Ral 5005 Signal Blue
--	----------------	---	---------------------	---	----------------------	---	--------------------	---	------------------------	---	----------------------

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.