

## Recirculation Fume Cupboard

Kinetic recirculation fume cupboards are available in standard widths of 1200mm, 1500mm, 1800mm and 2000mm.

Featuring aerodynamic zintec mild steel outer carcass with epoxy polyester powder coated finish, 316 grade stainless steel front sill, Trespa chambers and baffles, two independent toughed glass sliding sashes including fail safe mechanism, fluorescent light fittings and ventilated storage cupboard ideal for storing acids, solvents or flammables.

Fitted with pre-filters to remove particulates, carbon filters to absorb chemical vapours and centrifugal fan units to re-circulate air back into the laboratory.

Kinetic recirculation fume cupboards have programmable control panels that meet EN14175 incorporating sash high alarms, airflow alarms and fan stop/start.

Kinetic also offer bespoke external and internal widths, depth and heights and a varied quantity and type of filters dependent on the size and use of fume cupboard.



Re-circulation fume cupboard



Re-circulation fume cupboard



Re-circulation fume cupboard

## Colour Range

### Steelwork



Ral 9010  
White



Ral 7035  
Light Grey

### Stainless Steel



Stainless  
Steel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.