

## Semi-Mobile Ducted Fume Cupboard

Kinetic semi-mobile ducted fume cupboards are compliant with CLEAPSS Guide G9 – Fume Cupboards in Schools (Revision of DfEE Building Bulletin 88), BS EN 14175 and BREEAM (ENE11 & HEA17)

Kinetic semi-mobile ducted fume cupboards are designed to protect students and teachers from the fumes associated with GCSE and A Level curriculum activities and all cupboards are CLEAPSS tested.

Kinetic semi-mobile ducted fume cupboard is movable by up to 2 metres, locking into position and offers 360 degree visibility which is ideal for demonstrations,

Kinetic semi-mobile ducted fume cupboards are available in standard widths of 1000mm and 1200 mm. Featuring Trespa banded work surface, toughened glass sliding sash, built in storage compartment, integral high level lighting and airflow control panel.

### Turnkey Installation

Kinetic can also offer full turnkey installation packages on fume cupboards including a full PVC extraction system and fan unit.



Semi-mobile ducted fume cupboard

## Colour Range

### Steelwork



Ral 9010  
White



Ral 7035  
Light Grey

### Stainless Steel



Stainless  
Steel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.