



## Docking Station

Kinetic docking stations are a lockable enclosure that houses gas, power, drainage and water connections which allows quick release when the mobile fume cupboard is not in use or is being moved to a different location.

Kinetic docking stations can be mounted onto walls, panels and cupboard units, and feature quick connection fittings for water and gas, a Vulcathene drain adaptor and a 240V 13A single socket outlet.



Mobile fume  
cupboard docking  
station

## Colour Range

### Steelwork



Ral 9010  
White



Ral 7035  
Light Grey

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.