



Laboratory Suspended Unit

Kinetic laboratory suspended units can be installed onto the appropriate laboratory steelwork with the ability of being removed if purpose alters or for access to rear services or cleaning. The laboratory furniture range is available in a bespoke choice of sizes and colours.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit

Associated Products



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop

Dimensions

Overall Height	500 - 800mm
Length	400 - 1200mm
Carcass Depth	400 - 600mm

Hardware















Hinges	165° opening hinge
Drawer Sides	White epoxy coated metal drawer sides

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDf)

 White K101PE	 Grey U112PE
---	---

Laboratory Fronts - Laminate Faced MDF (LFMDf)

 Aluminium D1309WF	 Stone Oak D2217PR	 Natural Maple D375PR	 Beech D381PR	 White K101PE	 Grey U112PE
 Blue U125PE	 Yellow U131PE	 Orange U136PE	 Coral Red U148PE	 Anthracite U164PE	 Black U190PE
 Lime U2504PE	 Caribbean Blue U2510PE	 Magnolia U520PE			

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.