



Laboratory Suspended Unit

Kinetic laboratory suspended units can be installed onto the appropriate laboratory steelwork with the ability of being removed if purpose alters or for access to rear services or cleaning. The laboratory furniture range is available in a bespoke choice of sizes and colours.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit

Associated Products



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop

Dimensions

Overall Height	500 - 800mm
Length	400 - 1200mm
Carcass Depth	400 - 600mm

Hardware

Hinges	165° opening hinge
Drawer Sides	White epoxy coated metal drawer sides

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDf)

	White K101PE		Grey U112PE
--	------------------------	---	-----------------------

Laboratory Fronts - Laminate Faced MDF (LFMDf)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.