



## Laboratory Suspended Unit

Kinetic laboratory suspended units can be installed onto the appropriate laboratory steelwork with the ability of being removed if purpose alters or for access to rear services or cleaning. The laboratory furniture range is available in a bespoke choice of sizes and colours.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit

## Associated Products



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in A frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop



Suspended units in C frame support system with Trespa Toplab worktop

## Dimensions

Overall Height	500 - 800mm
Length	400 - 1200mm
Carcass Depth	400 - 600mm

## Hardware














Hinges	165° opening hinge
Drawer Sides	White epoxy coated metal drawer sides

## Colour Range

### Laboratory Carcass - Melamine Faced MDF (MFMDf)

 <b>White</b> K101PE	 <b>Grey</b> U112PE
---	---

### Laboratory Fronts - Laminate Faced MDF (LFMDf)

 <b>Aluminium</b> D1309WF	 <b>Stone Oak</b> D2217PR	 <b>Natural Maple</b> D375PR	 <b>Beech</b> D381PR	 <b>White</b> K101PE	 <b>Grey</b> U112PE
 <b>Blue</b> U125PE	 <b>Yellow</b> U131PE	 <b>Orange</b> U136PE	 <b>Coral Red</b> U148PE	 <b>Anthracite</b> U164PE	 <b>Black</b> U190PE
 <b>Lime</b> U2504PE	 <b>Caribbean Blue</b> U2510PE	 <b>Magnolia</b> U520PE			

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.