



## Laboratory Mobile Unit

Kinetic laboratory mobile units are fixed on swivel castors allowing them to be transported from one area to another. The laboratory furniture range is available in a bespoke choice of sizes and colours and can accommodate Grathnell trays.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit




Triple tray unit with doors



Triple open tray unit



Double tray unit with doors



Double open tray unit

## Associated Products



Mobile units in A frame support system with Trespa Toplab worktop



Mobile units in A frame support system with Trespa Toplab worktop



Mobile units in C frame support system with Trespa Toplab worktop



Mobile units in C frame support system with Trespa Toplab worktop



Mobile units in T frame support system with Trespa Toplab worktop



Mobile units in T frame support system with Trespa Toplab worktop

## Dimensions

Overall Height	500 - 827mm
Length	400 - 1000mm
Carcass Depth	400 - 600mm

## Hardware

Hinges	165° opening hinge 270° exposed axle hinge (Tray)
Feet	Swivel castors
Drawer Sides	White epoxy coated metal drawer sides
Trays	Gratnell trays in optional depth and colour

## Colour Range


### Laboratory Carcass - Melamine Faced MDF (MFMDF)

	<b>White</b> K101PE		<b>Grey</b> U112PE
----------------------------------------------------------------------------------	------------------------	-----------------------------------------------------------------------------------	-----------------------

### Laboratory Fronts - Laminate Faced MDF (LFMDF)

	<b>Aluminium</b> D1309WF		<b>Stone Oak</b> D2217PR		<b>Natural Maple</b> D375PR		<b>Beech</b> D381PR		<b>White</b> K101PE		<b>Grey</b> U112PE
	<b>Blue</b> U125PE		<b>Yellow</b> U131PE		<b>Orange</b> U136PE		<b>Coral Red</b> U148PE		<b>Anthracite</b> U164PE		<b>Black</b> U190PE
	<b>Lime</b> U2504PE		<b>Caribbean Blue</b> U2510PE		<b>Magnolia</b> U520PE						

### Gratnell Trays

	<b>Charcoal Grey</b>		<b>Light Grey</b>		<b>Black</b>		<b>Yellow</b>		<b>Red</b>		<b>Blue</b>
	<b>Green</b>		<b>Maroon</b>								

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.