



Laboratory Mobile Unit

Kinetic laboratory mobile units are fixed on swivel castors allowing them to be transported from one area to another. The laboratory furniture range is available in a bespoke choice of sizes and colours and can accommodate Grathnell trays.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit



Triple tray unit with doors



Triple open tray unit



Double tray unit with doors



Double open tray unit

Associated Products



Mobile units in A frame support system with Trespa Toplab worktop



Mobile units in A frame support system with Trespa Toplab worktop



Mobile units in C frame support system with Trespa Toplab worktop



Mobile units in C frame support system with Trespa Toplab worktop



Mobile units in T frame support system with Trespa Toplab worktop



Mobile units in T frame support system with Trespa Toplab worktop

Dimensions


Overall Height	500 - 827mm
Length	400 - 1000mm
Carcass Depth	400 - 600mm

Hardware





Hinges	165° opening hinge 270° exposed axle hinge (Tray)
Feet	Swivel castors
Drawer Sides	White epoxy coated metal drawer sides
Trays	Gratnell trays in optional depth and colour

Colour Range




Laboratory Carcass - Melamine Faced MDF (MFMDF)

 White K101PE	 Grey U112PE
---	---

Laboratory Fronts - Laminate Faced MDF (LFMDF)

 Aluminium D1309WF	 Stone Oak D2217PR	 Natural Maple D375PR	 Beech D381PR	 White K101PE	 Grey U112PE
 Blue U125PE	 Yellow U131PE	 Orange U136PE	 Coral Red U148PE	 Anthracite U164PE	 Black U190PE
 Lime U2504PE	 Caribbean Blue U2510PE	 Magnolia U520PE			

Gratnell Trays

 Charcoal Grey	 Light Grey	 Black	 Yellow	 Red	 Blue
 Green	 Maroon				

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.