



## Laboratory Wall Unit


Kinetic laboratory wall units are fixed onto walls, usually above laboratory workspaces to provide additional laboratory storage. The laboratory furniture range is available in a bespoke choice of sizes, colours and door types.



Double 2 door wall unit with solid doors



Single 1 door wall unit with solid door



Double 2 door wall unit with glazed doors



Single 1 door wall unit with glazed door



Double 2 door wall unit with glass doors



Single 1 door wall unit with glass door



Double 2 door wall unit with glass sliding doors

## Dimensions

Overall Height	400 - 600mm
Length	400 - 1000 mm
Carcass Depth	300 - 400 mm

## Hardware

Hinges	165° opening hinge Glass door hinge Glass sliding door track
--------	--

## Colour Range

### Laboratory Carcass - Melamine Faced MDF (MFMDF)

	<b>White</b> K101PE		<b>Grey</b> U112PE
--	------------------------	---	-----------------------

### Laboratory Fronts - Laminate Faced MDF (LFMDF)

	<b>Aluminium</b> D1309WF		<b>Stone Oak</b> D2217PR		<b>Natural Maple</b> D375PR		<b>Beech</b> D381PR		<b>White</b> K101PE		<b>Grey</b> U112PE
	<b>Blue</b> U125PE		<b>Yellow</b> U131PE		<b>Orange</b> U136PE		<b>Coral Red</b> U148PE		<b>Anthracite</b> U164PE		<b>Black</b> U190PE
	<b>Lime</b> U2504PE		<b>Caribbean Blue</b> U2510PE		<b>Magnolia</b> U520PE						

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.