

## Laboratory Grade Trespa TopLab®PLUS Worktop

TopLab®PLUS is a specified material that is chemically resistant, durable and easy to clean. Kinetic manufacture laboratory worktops for a number of different laboratory settings. The range is available in a bespoke choice of sizes and colours.



Trespa worktop with top 2x2mm chamfer



Trespa Toplab plus worktop



Trespa worktop with top 2x2mm chamfer



Trespa worktop with top 2x2mm chamfer



Trespa worktop with top & bottom 2x2mm chamfer



Trespa worktop with top & bottom 2x2mm chamfer



Trespa worktop with top & bottom 2x2mm chamfer



Trespa worktop with square edge



Trespa worktop with square edge



Trespa worktop with square edge



Trespa worktop with 7mm radiused edge



Trespa worktop with 7mm radiused edge



Trespa worktop with 7mm radiused edge

## Dimensions

Thickness	16 or 20mm
Depth	600 - 900mm

## Colour Range

### Laboratory Worktop - Trespa Toplab® Plus



White



Pastel Grey



Silver Grey



Slate Grey



Glacier Blue



Pastel Grey /  
Silver



Black

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.