



Velstone Laboratory Worktop

Velstone is a solid surface that is fire retardant, chemical resistant and easy to clean. Kinetic manufacture laboratory worktops for a number of different laboratory and educational settings. The range is available in a bespoke choice of sizes and colours.



Velstone worktop with pencil rounded edge



Velstone worktop



Velstone worktop with pencil rounded edge



Velstone worktop with pencil rounded edge

Dimensions

Thickness	25mm
Depth	600, 750, 900mm

Colour Range

Laboratory Worktop - Velstone



Soft White



Silver Ash



Packice



Silver Dust

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.