




## Underslung Sink Bowls

Kinetic offer a range of underslung epoxy sink bowls that are designed to be easily fitted to all types of worktop. The range is available in a choice of sizes, shapes and colours.




Underslung bowl with fluted drainer on sink unit




Underslung bowl dimensions




Underslung bowl without fluted drainer on sink unit




Underslung bowl with fluted drainer on sink unit



Underslung bowl without fluted drainer on sink unit



Underslung bowl with fluted drainer on sink unit




Underslung bowl without fluted drainer on sink unit



Underslung bowl with fluted drainer on sink unit




Underslung bowl without fluted drainer on sink unit




Underslung bowl mounting detail




Underslung bowl mounting detail




Underslung bowl in worktop




Underslung bowl in worktop




Underslung bowl in worktop



Underslung bowl in worktop



Underslung bowl in worktop



Underslung bowl in worktop

## Sizes

Code	Length	Width	Depth
	mm	mm	mm
SES2	160	160	150
SVS3	250	200	150
SDS5	300	200	185
SA5	356	254	152
SQ6	300	300	200
SS7	350	300	150
SS6	350	300	200
SKL7	380	350	200
SN9	400	300	200
SM8 (W)	450	368	130
SM8	450	356	200
SM8 (X)	450	356	280
SR12	550	365	250
SGL10	400	400	300
SGR11	500	400	300
SEL16	600	450	295
SPL14	610	406	203
SVL15	635	381	254

## Colour Range

### Sink - Epoxy Resin



Grey



White



Black

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.