

Laboratory Water Outlets

Kinetic offer a range of mountable laboratory water outlets suitable for any laboratory environment.

The laboratory outlets can be mounted on laboratory benches, laboratory worktops, walls, panels, overhead units and under shelves.

Available with colour coded handles and marked as per the EN13792 standard.

The outlet is white powder coated as standard but alternative colours are available.

Available as

Water potable, cold (WPC)

Water potable, hot (WPH)

Water non-potable, cold (WNC)

Water non-potable, hot (WHH)



1 way worktop mounted



1 way reagent mounted



1 way reagent or wall mounted



1 way wall mounted



Colour Range

Other Colour Options



Grey Powder
Coated



Polished
Chrome



Satin Chrome



Satin Nickel

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.