

## Laboratory Mixed Water Taps

Kinetic offer a range of mountable mixed water taps suitable for any laboratory environment.

The taps can be mounted on laboratory benches and are available in swanneck or pillar angle fitting.

Available in lever action or as a standard handle which are colour coded and marked as per the EN13792 standard.

The outlet is white powder coated as standard but alternative colours are available.


Available as

Water potable, cold (WPC)

Water potable, hot (WPH)

Water non-potable, cold (WNC)

Water non-potable, hot (WHH)





Hot & cold mixer tap with hand wheels & removable nozzle



Hot & cold mixer tap with wrist action levers & removable nozzle



Hot & cold monobloc tap with hand wheels & removable nozzle



Hot & cold monobloc tap with wrist action levers & removable nozzle



Wall mounted hot & cold mixer tap with hand wheels & removable nozzle



Wall mounted hot & cold mixer tap with wrist action levers & removable nozzle



Hot & cold monobloc tap with single handed wrist action lever & removable nozzle



Hot or cold twin inline tap with hand wheels



Hot or cold twin tap with hand wheels



Hot or cold triple tap with hand wheels

# Colour Range

## Other Colour Options



Grey Powder Coated



Polished Chrome



Satin Chrome



Satin Nickel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.