

## Laboratory Precise Flow Technical Gas Needle Outlets

Kinetic offer a range of laboratory technical gas needle outlets that are suitable for any laboratory environment.

The outlets are available as single or twin and can be mounted on laboratory benches, laboratory worktops, walls, panels, overhead units and reagent shelves.

Available with colour coded handles and marked as per the EN13792 standard.

The outlet is white powder coated as standard but alternative colours are available.

Available as

Hydrogen (H<sub>2</sub>)

Compressed Air (CA)

Oxygen (O<sub>2</sub>)

Nitrogen (N<sub>2</sub>)

Carbon Dioxide (CO<sub>2</sub>)

Argon (Ar)

Hellium (He)

Dinitrogen monoxide (N<sub>2</sub>O)



1 way worktop mounted



1 way reagent mounted



1 way reagent or wall mounted



2 way wall mounted



1 way reagent mounted



1 way wall mounted

## Colour Range

### Other Colour Options



Grey Powder  
Coated



Polished  
Chrome



Satin Chrome



Satin Nickel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.