

## Laboratory Push Turn Burning Gas Outlets

Kinetic offer a range of laboratory burning gas needle outlets that are suitable for any laboratory environment.

The outlets are available as single or twin and can be mounted on laboratory benches, laboratory worktops, walls, panels, overhead units and reagent shelves.

Available with colour coded handles and marked as per the EN13792 standard.

The outlet is white powder coated as standard but alternative colours are available.

Available as

Natural Gas (G)

Liquified Petrol Gas (LPG)

Butane (C<sub>4</sub>H<sub>10</sub>)

Propane (C<sub>3</sub>H<sub>8</sub>)

Acetylene (C<sub>2</sub>H<sub>2</sub>)

Methane (CH<sub>4</sub>)



1 way worktop mounted



1 way reagent mounted



1 way reagent or wall mounted



2 way wall mounted



1 way wall or wall mounted



1 way wall mounted

# Colour Range

## Other Colour Options



Grey Powder Coated



Polished Chrome



Satin Chrome



Satin Nickel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.