

## Laboratory Vacuum Outlets

Kinetic offer a range of laboratory vacuum outlets that are suitable for special laboratory processes.

The outlets are available as single or twin and can be mounted on Laboratory benches, Laboratory worktops, walls, panels, overhead units and reagent shelves.

Available with colour coded handles and marked as per the EN13792 standard.

The outlet is white powder coated as standard but alternative colours are available.

Available as

Low Vacuum – 100kPa to 0,1 kPa (V)

Fine Vacuum – 0,1kPa to 0,001 kPa (VF)


High Vacuum – 0,001kPa to 0,0000001 kPa (VH)



1 way worktop mounted



2 way inline reagent mounted



1 way reagent or wall mounted



1 way wall mounted

## Colour Range

### Other Colour Options



Grey Powder  
Coated



Polished  
Chrome



Satin Chrome



Satin Nickel

Colours on this page may not be a true representation.  
For accurate colour samples, please contact us.