



Mobile Unit

Kinetic mobile units are fixed on swivel castors allowing them to be transported from one area to another. The range is available in a bespoke choice of sizes and colours and can accommodate Gratnell trays.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit



Triple tray unit with doors



Triple open tray unit

Dimensions

Overall Height	500 - 827mm
Length	400 - 1000mm
Carcass Depth	400 - 600mm

Hardware

Hinges	165° opening hinge 270° exposed axle hinge (Tray) Glass door hinge Glass sliding door track
Feet	Swivel castors
Drawer Sides	White epoxy coated metal drawer sides
Trays	Gratnell trays in optional depth and colour

Colour Range


Carcass - Melamine Faced MDF (MFMDf)

	White K101PE		Grey U112PE		Beech D381PR
--	------------------------	---	-----------------------	---	------------------------

Laboratory Fronts - Laminate Faced MDF (LFMDf)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Gratnell Trays

	Charcoal Grey		Light Grey		Black		Yellow		Red		Blue
	Green		Maroon								

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.