



Pedestal Unit

Kinetic pedestal units stand on adjustable feet allowing the unit to be positioned and levelled to support the worktop. The range is available in a bespoke choice of sizes, colours and can accommodate Gratnell trays.



Double 2 door unit



Double 2 door 1 drawer unit



Double 2 door 2 drawer unit



Double 4 drawer unit



Single 1 door unit



Single 1 door 1 drawer unit



Single 4 drawer unit



Triple tray unit with doors



Triple open tray unit



Double sink unit



Single sink unit

Dimensions

Overall Height	600 - 1000mm
Length	400 - 1200mm
Carcass Depth	400 - 600mm

Hardware

Hinges	165° opening hinge 270° exposed axle hinge (Tray)
Feet	Height adjustable plinth feet
Drawer Sides	White epoxy coated metal drawer sides
Trays	Gratnell trays in optional depth and colour

Colour Range



Carcass - Melamine Faced MDF (MFMDf)

	White K101PE		Grey U112PE		Beech D381PR
--	------------------------	---	-----------------------	---	------------------------

Laboratory Fronts - Laminate Faced MDF (LFMDf)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Gratnell Trays

	Charcoal Grey		Light Grey		Black		Yellow		Red		Blue
	Green		Maroon								

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.