

Pedestal

Kinetic Pedestal support system is fixed into position and directly supports the worktop.


The under bench unit stands on a plinth in a fixed position and can vary in length.

The storage spaces can be customised using a combination of drawers, shelves and cupboard spaces.

The Pedestal support system is available in a range of colours.



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop



Pedestal units with Beech laminate worktop

Available On

- Melamine Faced MDF
- Laminate MDF

Colour Range

Carcass - Melamine Faced MDF (MF MDF)



Laboratory Fronts - Laminate Faced MDF (LF MDF)



Laboratory Worktop - Laminate



Colours on this page may not be a true representation.
For accurate colour samples, please contact us.