



A Frame

Kinetic A Frame support system allows independent worktop support without the use of under bench storage.

The span frames and the low level back rails can vary in length and will form a series of individual spaces. Allowing the client the option of adding suspended units, free standing units or mobile units.

The storage spaces can also be customised to take clients own equipment such as fridges and freezers or can be left open to provide knee spaces for write-up areas.

Steel frame are epoxy powder coated and are available in a range of colours.



A Frame support system with Beech laminate worktop




A Frame support system with Beech laminate worktop



A Frame support system with Beech laminate worktop



A Frame support system with Beech laminate worktop



A Frame support system with Beech laminate worktop



A Frame support system with Beech laminate worktop



A Frame support system with Beech laminate worktop

Colour Range

Laboratory Worktop - Laminate



Aluminium
D1309WF



Stone Oak
D2217PR



Natural Maple
D375PR



Beech
D381PR



White
K101PE



Grey
U112PE



Blue
U125PE



Yellow
U131PE



Orange
U136PE



Coral Red
U148PE



Anthracite
U164PE



Black
U190PE



Lime
U2504PE



Caribbean Blue
U2510PE



Magnolia
U520PE

Steelwork



Ral 9010
White



Ral 7035
Light Grey



Ral 7012
Basalt Grey



Ral 7021
Dark Grey



Ral 5003
Sapphire Blue

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.