



Tall Storage Unit

Kinetic tall storage units are a self-supporting unit supported on a fixed plinth with encase adjustable feet for levelling. The range is available in a bespoke choice of sizes and colours.



Double 2 door tall storage unit



Double open tall storage unit



Double 2 door 4 drawer tall storage unit



Double 4 drawer open tall storage unit



Single 1 door tall storage unit



Single open tall storage unit



Single 4 drawer open tall storage unit



Single 1 door 4 drawer tall storage unit

Dimensions

Overall Height	1500 - 1750mm
Length	400 - 1000mm
Carcass Depth	400 - 600mm

Hardware














Hinges	165° opening hinge
Feet	Encased height adjustable feet
Drawer Sides	White epoxy coated metal drawer sides

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDF)

 White K101PE	 Grey U112PE
---------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------

Laboratory Fronts - Laminate Faced MDF (LFMDF)

 Aluminium D1309WF	 Stone Oak D2217PR	 Natural Maple D375PR	 Beech D381PR	 White K101PE	 Grey U112PE
 Blue U125PE	 Yellow U131PE	 Orange U136PE	 Coral Red U148PE	 Anthracite U164PE	 Black U190PE
 Lime U2504PE	 Caribbean Blue U2510PE	 Magnolia U520PE			

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.