



Tall Storage Unit

Kinetic tall storage units are a self-supporting unit supported on a fixed plinth with encase adjustable feet for levelling. The range is available in a bespoke choice of sizes and colours.



Double 2 door tall storage unit



Double open tall storage unit



Double 2 door 4 drawer tall storage unit



Double 4 drawer open tall storage unit



Single 1 door tall storage unit



Single open tall storage unit



Single 4 drawer open tall storage unit



Single 1 door 4 drawer tall storage unit

Dimensions

Overall Height	1500 - 1750mm
Length	400 - 1000mm
Carcass Depth	400 - 600mm

Hardware
















Hinges	165° opening hinge
Feet	Encased height adjustable feet
Drawer Sides	White epoxy coated metal drawer sides

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDF)

	White K101PE		Grey U112PE
--	------------------------	---	-----------------------

Laboratory Fronts - Laminate Faced MDF (LFMDF)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.