



Wall Unit

Kinetic wall units are fixed onto walls, usually above workspaces to provide additional storage. The range is available in a bespoke choice of sizes, colours and door types.



Double 2 door wall unit with solid doors



Single 1 door wall unit with solid door



Double 2 door wall unit with glazed doors



Single 1 door wall unit with glazed door



Double 2 door wall unit with glass doors



Single 1 door wall unit with glass door



Double 2 door wall unit with glass sliding doors

Dimensions

Overall Height	400 - 600mm
Length	400 - 1000 mm
Carcass Depth	300 - 400 mm

Hardware
















Hinges	165° opening hinge Glass door hinge Glass sliding door track
--------	--

Colour Range

Laboratory Carcass - Melamine Faced MDF (MFMDF)

	White K101PE		Grey U112PE
--	------------------------	---	-----------------------

Laboratory Fronts - Laminate Faced MDF (LFMDF)

	Aluminium D1309WF		Stone Oak D2217PR		Natural Maple D375PR		Beech D381PR		White K101PE		Grey U112PE
	Blue U125PE		Yellow U131PE		Orange U136PE		Coral Red U148PE		Anthracite U164PE		Black U190PE
	Lime U2504PE		Caribbean Blue U2510PE		Magnolia U520PE						

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.