




Sit on Sink Bowls

Kinetic offer a range of sit on epoxy sink bowls that are designed to be easily fitted to all types of worktop. The range is available in a choice of sizes, shapes and colours.




Sit on bowl on sink unit



Sit on bowl dimensions



Sit on bowl on sink unit



Sit on bowl on sink unit



Sit on bowl on sink unit




Sit on bowl mounting detail



Sit on bowl mounting detail



Sit on bowl mounting detail




Sit on bowl in worktop



Sit on bowl in worktop



Sit on bowl in worktop



Sit on bowl in worktop



Sit on bowl in worktop



Sit on bowl in worktop

Sizes

Code	Length	Width	Depth
	mm	mm	mm
SGSS-DO	256	108	150
SDS5-DO	300	200	200
SS6-DO	350	300	200
SM8-DO	450	365	200
SGL10-DO	400	400	300
SEL16-DO	600	450	300

Colour Range

Laboratory Worktop - Epoxy Resin



Grey



White



Black

Colours on this page may not be a true representation.
For accurate colour samples, please contact us.