

General fixtures/ furnishings/ equipment

N10 General fixtures/ furnishings/ equipment

MODULAR TABLES

- Manufacturer: Kinetic Laboratories Ltd, Unit 7-9 GB Business Park, Wiltshire Road Hull, HU4 6PA
 - Product reference: Modular Free Standing table frame.
- Height (overall): 900mm
- Width (overall): 750mm
- Length (overall): 1500mm
- Worktop: Trespa Toplab Base, Colour TBC
- Adjustable feet: 20mm adjustment, max load capacity 240kg.
- Metalwork materials: 40x40x16swg ERW steel tube, 40x20x16swg ERW steel tube
- Ironmongery/ Accessories: Slot plates
- Other requirements: None.

INSTALLATION

- Programming: Do not install tables before building is weathertight, wet trades have finished their work, wall and floor finished are complete, and the building is well dried out.
- Modifications: Do not cut, plane or sand prefinished components except where shown on drawings.
- Fixing: Secure components using methods and fasteners recommended by the manufacturer.